	8h	Nursing/O ther	12h	14h	Nursing/ Other	18h	22h	Nursing/o ther	2 h	6 h
VENTILATOR 1										
The patient is on a life support ventilator										
The ventilator is plugged and it is confirmed on the screen										
The ventilator continues on the battery if unplugged										
Disconnection alarm rings in <20s										
The ventilator is not in its transport bag										
VENTILATOR 2										
A second ventilator »life support » is in the room										
The 2nd ventilator is plugged and it is confirmed on the screen.										
The tubes are connected to the 2 nd ventilator										
The ventilator starts ti ventilate if switched on										
MANUAL RESUSCITATOR										
An ambu Bag is in the room, ready to be used										
CANULA										
A 2 nd canula and a connexion hose are available in the room										
NURSE NAME										